<u>User Control Panel (0 new messages)</u> • <u>View your posts</u> • <u>Logout [potto216]</u>

Skip to content

- FAO
- Members
- Search
- Board index \(\cdot\) Platforms \(\cdot\) Sonix RP
- Change font size E-mail friend Print view

Post a reply

Search

First unread post • 3 posts • Page 1 of 1

Correlating Frame Rate Values from Calculations

- Edit post
- Report this post
- Reply with quote

Correlating Frame Rate Values from Calculations

by **potto216** » Sun Nov 13, 2011 5:32 pm

A question about the frame rate computation in Sonix RP. When we setup the collect for a line density of 64 and a sector size of 50% with a depth of 2.5cm we get an output file with a width of 33. If I compute the theoretical frame rate I get:

Time_per_line=2*(2.5 cm/1540 m/s)=32.46 usFrame_rate=1/(33*32.46 us)=933 fps

But the actual frame rate is 532fps. Are my calculations incorrect or is there a processing delay between each from collect? If so we need to know if the delay occurs at the end of the frame collect or if it is inter-spaced between each collected line?

Thanks!

potto216

Posts: 21

Joined: Thu Feb 18, 2010 4:52 am

• Private message

Top

- Report this post
- Reply with quote

Re: Theoretical frame rate does not match Sonix RP frame rat

by <u>reza</u> » Tue Apr 10, 2012 10:05 pm

This is expected, the imaging frame rate will always be smaller than the theoritical frame rate. This is because there is always a minimum switching delay in the hardware when you switch from one scan line to the next. This switching can not happen instantly. This delay increases your line duration which results in lower frame rates.

Actual line duration can be calculated by dividing your frame rate by the number of scan lines which in your case is

line duration = 1/532/33

Hope this helps

Reza Zahiri Azar, PhD | Research Manager Ultrasonix Medical Corporation | #130 - 4311 Viking Way, BC, Canada V6V2K9 T 604.279.8550 x 185 | Skype reza.zahiri http://www.ultrasonix.com

<u>reza</u>

Posts: 362

Joined: Tue Dec 13, 2005 9:34 pm **Location:** University of British Columbia

- Private message
- Website
- AIM

Top

- Report this post
- Reply with quote

Re: Theoretical frame rate does not match Sonix RP frame rat

by **kris** » Wed Jun 20, 2012 9:56 pm

Note the probes.xml has a parameter called minlineduration, this ensures that the PRF never gets too high in B mode imaging (safety).

For example, it's 60us for L14-5/38, which explains your \sim 530 fps.

Kris Dickie | Managing Director Ultrasonix Medical Corporation | #130 - 4311 Viking Way, BC, Canada V6V2K9 T 604.279.8550 x 104 | Skype krisdickie | www.ultrasonix.com

<u>kris</u>

Site Admin

Posts: 731

Joined: Mon Dec 12, 2005 11:02 pm **Location:** Ultrasonix Medical Corporation

- Private message
- Website

	\mathbf{O}	n
J	U	ν

Display posts from previous: All posts ▼ Sort by Post time ▼ Ascending ▼ Go									
	Display posts from previous:	All posts	•	Sort by	Post time	▼	Ascending	•	Go

Post a reply

3 posts • Page 1 of 1

Return to Sonix RP

Jump to:	Sonix RP	▼ Go
1		

Who is online

Users browsing this forum: potto216 and 0 guests

• <u>Subscribe topic</u> • <u>Bookmark topic</u> • <u>Bump topic</u> • <u>The team</u> • <u>Delete all board cookies</u> • All times are UTC

Powered by phpBB® Forum Software © phpBB Group Cerulean Style by Daniel St. Jules of Gamexe.net